

PATIENT INFORMATION SHEET FOR THE COLLECTION OF FAECES (STOOL) M/C/S

METHOD OF COLLECTION:

1. You are given one brown lid, white plastic jar/container.
2. Pass the faeces into a clean plastic container and using the scoop located on the lid of the brown top jar/ container, transfer a small portion of faeces into a brown-top jar/container.
3. Place the lid on the jar/container and tighten the lid firmly to prevent leakage. Discard the remainder of the sample into the toilet.
4. Label the jar/container with your name, date of birth, date and time of collection. Deliver to the nearest Medlab Collection Centre or to your Doctor's surgery for transport arrangements.

For further inquiries or additional information regarding this testing procedure you can either speak directly with the attending pathology collector or by contacting the main laboratory Customer Service Department on 1300 633 522.

Issue No / Document ID:	2.0.1 / 25220	Issue Date:	2/12/2020
Prepared by:	Maja Trajkova	Authorised by:	Dr K.P. Singh